





















Hex Tool Application

Hex Tool Used		Compatible Parts Used With Specified Hex Tool					
1.25 Hex Tool		All Cover Screws	All Healing Collars	All Ball Abutments	All Screw-Receiving Abutments	All Abutment Screws	
							
2.0 Hex Tool		ScrewDirect 3030 & 3037 Implants		Legacy3 3.2mm implants 3.0 platform (8532)	All ScrewRedirect Implants		
							
2.5 Hex Tool		All FixtureMounts	Screwplant 1037 & 47 Screwplus 2030 & 37	ScrewDirect 3047 & 3057 Implants	Legacy(3) 3.5 & 4.5 platforms	All ScrewIndirect Implants	
							
3.0 Hex Tool		ScrewPlant/Plus Implants 1057, 2047 & 2057		Legacy(3) 5.7 platform			
