



Enhanced Esthetics

Introducing ease and economy into your esthetic zone restorations



Zirconia Abutments

Range of options

- three angulations (0, 8 or 15°)
- two collar heights (1 or 2mm)
- four connections (Spectra-System; industry-standard internal hex, tri-lobe or octagon)

All-in-One Packaging

- abutment
- system-specific fixation screw

Contoured Zirconia Head

- contoured margins and collar height choices for minimal preparation
- thickness allows for substantial customization without reducing core strength

Titanium Base

- provides structural support for zirconia
- creates a titanium screw seat for reduced possibility of fracture during abutment attachment
- engages implant for a secure titanium-to-titanium connection
- gold anodized for simulation of natural tooth color

System-Specific Fixation Screw

- all compatible with standard 1.25mmD (.050in) hex tool
- RePlant® compatible with Nobel Biocare's Uni-Grip™ Driver
- Swish™ compatible with Straumann's SCS Screwdriver¹

¹Except for 3.7mmD Swish™ abutments

0°, 8° & 15° Zirconia Contoured Abutments

Includes Dual-Grip™ Fixation Screw

Code	Platform	Angulation	Collar Heights	Replacement Screw	Price
6035-31C	3.5mmD	0°	1mm	1000-23	\$120.00
6035-32C	3.5mmD	0°	2mm	1000-23	\$120.00
6035-51-08C	3.5mmD	8°	1mm	1000-23	\$120.00
6035-52-08C	3.5mmD	8°	2mm	1000-23	\$120.00
6035-51-15C	3.5mmD	15°	1mm	1000-23	\$120.00
6035-52-15C	3.5mmD	15°	2mm	1000-23	\$120.00
6035-30CB	3.5mmD	-	.5mm	1000-23	\$60.00

Code	Platform	Angulation	Collar Heights	Replacement Screw	Price
6043-31C	4.3mmD	0°	1mm	1000-22	\$120.00
6043-32C	4.3mmD	0°	2mm	1000-22	\$120.00
6043-51-08C	4.3mmD	8°	1mm	1000-22	\$120.00
6043-52-08C	4.3mmD	8°	2mm	1000-22	\$120.00
6043-51-15C	4.3mmD	15°	1mm	1000-22	\$120.00
6043-52-15C	4.3mmD	15°	2mm	1000-22	\$120.00
6043-30CB	4.3mmD	-	.5mm	1000-22	\$60.00

Code	Platform	Angulation	Collar Heights	Replacement Screw	Price
6050-31C	5.0mmD	0°	1mm	1000-22	\$120.00
6050-32C	5.0mmD	0°	2mm	1000-22	\$120.00
6050-51-08C	5.0mmD	8°	1mm	1000-22	\$120.00
6050-52-08C	5.0mmD	8°	2mm	1000-22	\$120.00
6050-51-15C	5.0mmD	15°	1mm	1000-22	\$120.00
6050-52-15C	5.0mmD	15°	2mm	1000-22	\$120.00
6050-30CB	5.0mmD	-	.5mm	1000-22	\$60.00



Straight Zirconia Contoured Abutments

Includes Fixation Screw [1000-01]

Code	Platform	Collar Heights	Price
1037-30CB	3.7mmD	.5mm	\$60.00
1037-31C	3.7mmD	1.5mm	\$120.00
1037-32C	3.7mmD	2.5mm	\$120.00
1047-30CB	4.7mmD	.5mm	\$60.00
1047-31C	4.7mmD	1.5mm	\$120.00
1047-32C	4.7mmD	2.5mm	\$120.00



ScrewPlant® -continued

8° & 15° Zirconia Contoured Abutments

Includes Fixation Screw [1000-01]

Code	Platform	Angulation	Collar Heights	Price
1037-51-08C	3.7mmD	8°	1.5mm	\$120.00
1037-51-15C	3.7mmD	15°	1.5mm	\$120.00
1037-52-08C	3.7mmD	8°	2.5mm	\$120.00
1037-52-15C	3.7mmD	15°	2.5mm	\$120.00
1047-51-08C	4.7mmD	8°	1.5mm	\$120.00
1047-51-15C	4.7mmD	15°	1.5mm	\$120.00
1047-52-08C	4.7mmD	8°	2.5mm	\$120.00
1047-52-15C	4.7mmD	15°	2.5mm	\$120.00



Legacy™

0°, 8° & 15° Zirconia Contoured Abutments

Price includes Fixation Screw with 1.25mmD Hex

Code	Platform	Angulation	Collar Heights	Replacement Screw	Price
8530-31C	3.0mmD	0°	1mm	1000-04	\$120.00
8530-32C	3.0mmD	0°	2mm	1000-04	\$120.00
8530-51-08C	3.0mmD	8°	1mm	1000-04	\$120.00
8530-52-08C	3.0mmD	8°	2mm	1000-04	\$120.00
8530-51-15C	3.0mmD	15°	1mm	1000-04	\$120.00
8530-52-15C	3.0mmD	15°	2mm	1000-04	\$120.00
8530-30CB	3.0mmD	-	.5mm	1000-04	\$60.00

Code	Platform	Angulation	Collar Heights	Replacement Screw	Price
8035-31C	3.5mmD	0°	1mm	1000-01	\$120.00
8035-32C	3.5mmD	0°	2mm	1000-01	\$120.00
8035-51-08C	3.5mmD	8°	1mm	1000-01	\$120.00
8035-52-08C	3.5mmD	8°	2mm	1000-01	\$120.00
8035-51-15C	3.5mmD	15°	1mm	1000-01	\$120.00
8035-52-15C	3.5mmD	15°	2mm	1000-01	\$120.00
8035-30CB	3.5mmD	-	.5mm	1000-01	\$60.00

Code	Platform	Angulation	Collar Heights	Replacement Screw	Price
8045-31C	4.5mmD	0°	1mm	1000-01	\$120.00
8045-32C	4.5mmD	0°	2mm	1000-01	\$120.00
8045-51-08C	4.5mmD	8°	1mm	1000-01	\$120.00
8045-52-08C	4.5mmD	8°	2mm	1000-01	\$120.00
8045-51-15C	4.5mmD	15°	1mm	1000-01	\$120.00
8045-52-15C	4.5mmD	15°	2mm	1000-01	\$120.00
8045-30CB	4.5mmD	-	.5mm	1000-01	\$60.00

